

GeneLab Standard Operating Procedure:  
Normalizing sequencing library pools without the use of  
iSeq sequencing output

*Sep 2021*

*Version 1.1*

## Document Revisions

Document Number	Revision Number	Date	Description of Changes
GL-SOP-6.4	1.0	Jan 2021	Original
GL-SOP-6.4	1.1	Sep 2021	Adjusted scope of procedure to incorporate non TruSeq libraries Changed procedure name to reflect that it can be used for various library kits Step – wording Added TapeStation, PicoGreen and Qubit consumables to “Reagents” Added missing steps to combine individual libraries

## Scope and Purpose

This SOP describes the steps taken by NASA GeneLab to normalize sequencing libraries in a pool.

## Equipment and Consumables

1. DNase- and RNase-free water
2. Ice
3. Ice bucket
4. low bind 1.5mL microtube
5. 96-well sterile plate with working capacity of 100-150u

## Reagents

1. TapeStation reagents (as described in GL-SOP-6.3)
2. Qubit reagents (as described in GL-SOP-4.1)
- Or
3. Quant-iT PicoGreen reagents (as described in GL-SOP-6.2)

## Procedure

1. Obtain average library fragment size from TapeStation D1000 following SOP GL-SOP-6.3. Check there is no and that the average size of libraries is within similar range  $\pm 50$  bp.
2. Convert each library concentration obtained from PicoGreen or Qubit measurement to molarity using average or individual fragment size in units of [bp].

Equation for converting dsDNA:

$$\text{conc nM} = \frac{(\text{conc in ng/uL}) \times 10^6}{(660 \text{ g/mol} \times \text{library size bp})}$$

3. Dilute each library to 20nM in DNase-/RNase-free water.
4. Dilute further down to required molarity and combine equal volumes of each of the library in to a 1.5mL low bind tube to create a library pool.
5. Dilute the library pool to loading concentration or store at  $-20^{\circ}\text{C}$  for iSeq and/or NovaSeq run (GL-SOP-7.1)

Sample #	Sample	Library PicoGreen conc. (ng/ul)	Libraries (nM)	20 nM	
				Sample Vol. to 20 nM in 20ul	Water ul
1	MGS_HLU IR_M9_D SKN_RNA_ALQ0	6.490	32.78	12.20	7.80

Figure 1: Example dilution of a sample.